

#### References

- 1 Wallentin, L., Hijazi, Z., Andersson, U., Alexander, J.H., De Caterina, R., Hanna, M. et al. (2014). Growth differentiation factor 15, a marker of oxidative stress and inflammation, for risk assessment in patients with atrial fibrillation: insights from the Apixaban for Reduction in Stroke and Other Thromboembolic Events in Atrial Fibrillation (ARISTOTLE) trial. *Circulation* 130, 1847–1858.
- 2 O'Brien, E.C., Simon, D.N., Thomas, L.E., Hylek, E.M., Gersh, B.J., Ansell, J.E. et al. (2015). The ORBIT bleeding score: a simple bedside score to assess bleeding risk in atrial fibrillation. *Eur Heart J* 36, 3258–3264.
- 3 Hijazi, Z., Oldgren, J., Lindbäck, J., Alexander, J.H., Connolly, S.J., Eikelboom, J.W. et al. (2016). The novel biomarker-based ABC (age, biomarkers, clinical history)-bleeding risk score for patients with atrial fibrillation: a derivation and validation study. *Lancet* 387, 2302–2311.
- 4 Wollert, K.C., Kempf, T., Peter, T., Olofsson, S., James, S., Johnston, N. et al. (2007). Prognostic value of growth-differentiation factor-15 in patients with non-ST-elevation acute coronary syndrome. *Circulation* 115, 962–971.
- 5 Wollert, K.C., Kempf, T., Lagerqvist, B., Lindahl, B., Olofsson, S., Allhoff, T. et al. (2007). Growth differentiation factor 15 for risk stratification and selection of an invasive treatment strategy in non ST-elevation acute coronary syndrome. *Circulation* 116, 1540–1548.
- 6 Damman, P., Kempf, T., Windhausen, F., van Straalen, J.P., Guba-Quint, A., Fischer, J. et al. (2014). Growth-differentiation factor 15 for long-term prognostication in patients with non-ST-elevation acute coronary syndrome: an Invasive versus Conservative Treatment in Unstable coronary Syndromes (ICTUS) substudy. *Int J Cardiol* 172, 356–363.
- 7 Bonaca, M.P., Morrow, D.A., Braunwald, E., Cannon, C.P., Jiang, S., Breher, S. et al. (2011). Growth differentiation factor-15 and risk of recurrent events in patients stabilized after acute coronary syndrome: observations from PROVE IT-TIMI 22. *Arterioscler Thromb Vasc Biol* 31, 203–210.
- 8 Wallentin, L., Lindhagen, L., Årnström, E., Husted, S., Janzon, M., Johnsen, S.P. et al. (2016). Early invasive versus non-invasive treatment in patients with non-ST-elevation acute coronary syndrome (FRISC-II): 15 year follow-up of a prospective, randomised, multicentre study. *Lancet* 388, 1903–1911.
- 9 Yancy, C.W., Jessup, M., Bozkurt, B., Butler, J., Casey, D.E., Drazner, M.H. et al. (2013). 2013 ACCF/AHA guideline for the management of heart failure: a report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines. *J Am Coll Cardiol* 62, e147–239.
- 10 Wollert, K.C., Kempf, T. (2012). Growth differentiation factor 15 in heart failure: an update. *Curr Heart Fail Rep* 9, 337–345.
- 11 Chan, M.M.Y., Santhanakrishnan, R., Chong, J.P.C., Chen, Z., Tai, B.C., Liew, O.W. et al. (2016). Growth differentiation factor 15 in heart failure with preserved vs. reduced ejection fraction. *Eur J Heart Fail* 18, 81–88.
- 12 Internal data (2016); data on file.
- 13 Elecsys® GDF-15 package insert (2017). Roche Diagnostics Documentation, Basel.

#### Published by

Roche Diagnostics International Ltd  
CH-6343 Rotkreuz  
Switzerland

© 2017

All trademarks mentioned enjoy legal protection.

[www.roche.com](http://www.roche.com)

## illuminating the hidden severity of cardiovascular disease with GDF-15

*A valuable biomarker to improve risk stratification and optimize patient management*



# Elecsys® GDF-15 Assay

## Elecsys® GDF-15 Assay – Intended Use

Immunoassay for the in vitro quantitative determination of Growth Differentiation Factor-15 (GDF-15) in human serum and plasma. The Elecsys GDF-15 assay is intended as an aid in risk stratification of patients with Acute Coronary Syndrome (ACS) or Chronic Heart Failure (CHF). The Elecsys GDF-15 assay is intended as an aid in risk prediction of major bleeding events of patients with Atrial Fibrillation (AF). The electrochemiluminescence immunoassay “ECLIA” is intended for use on Elecsys and **cobas e** immunoassay analyzers.<sup>13</sup>

## Elecsys® GDF-15 Assay – Value Propositions

### Predicting bleeding risk in atrial fibrillation patients with more confidence

- Bleeding risk is a leading concern when selecting oral anticoagulation (OAC) therapy for stroke prevention in patients with atrial fibrillation (AF)<sup>1,2</sup>
- GDF-15 is the strongest predictor in the ABC bleeding risk score.<sup>3</sup> It provides physicians with an improved understanding of each patient's bleeding risk profile and can inform treatment decisions<sup>3</sup>

### Better decisions in acute coronary syndrome


- GDF-15 is a strong independent prognostic marker to anticipate all-cause mortality and reinfarction in patients with acute coronary syndrome (ACS)<sup>4-7</sup>
- The addition of GDF-15 to TnT-hs can facilitate the selection of patients with non-ST elevation acute coronary syndrome (NSTEMI-ACS) who would benefit from early revascularization and more intensive medical therapies<sup>5,8</sup>

### Get the full picture in heart failure

- Heart failure (HF) is a complex syndrome, affecting not only the heart but many extra-cardiac tissues<sup>9,10</sup>
- GDF-15 is independently associated with mortality and non-fatal events in HF with reduced (HFrEF) or preserved ejection fraction (HFpEF)<sup>11</sup>
- Together with NT-proBNP and TnT-hs it provides a comprehensive clinical picture and helps to identify high-risk patients who benefit from intensified management<sup>11</sup>

## Elecsys® GDF-15 Assay – Technical features

### Lot to Lot Method Comparison according to CLSI EP09-A3




Lot to lot method comparison according to CLSI EP09-A3 between P02 and P03 lot on cobas e 411. Sample concentrations were between 450 and 19,944 pg/mL (P02 lot).<sup>12,13</sup>

Passing/Bablok  
 $y = 0.996x + 1.22$   
 $t = 0.976$

Linear regression  
 $y = 0.995x + 18.3$   
 $r = 1.00$

N = 138

### Elecsys® platform comparison according to CLSI EP09-A3



Platform method comparison according to CLSI EP09-A3 between cobas e 601 and cobas e 801. Sample concentrations were between 403 and 18,876 pg/mL (cobas e 601).<sup>12,13</sup>

Passing/Bablok  
 $y = 1.017x - 3.97$   
 $t = 0.989$


Linear regression  
 $y = 1.021x - 10.9$   
 $r = 1.00$

N = 143

## Proven Elecsys® assay performance – confidence in results

- Advanced ECL assay design and low inherent system variability
- High precision at across measure range

### Repeatability/Within-Run-Precision of Elecsys® GDF-15 immunoassay according to CLSI-EP05-A3



### Intermediate/Within-Lab Precision of Elecsys® GDF-15 immunoassay according to CLSI-EP05-A3




Fig. 3: Precision according to the EP5-A3 protocol of the Clinical and Laboratory Standards Institute (CLSI) on cobas e 411 and e 601.<sup>12,13</sup>

## Elecsys® GDF-15 test characteristics<sup>12,13</sup>

Testing time	18 min
Test principle	Sandwich principle
Calibration	2-point calibration
Traceability	The assay has been standardized by weighing recombinant GDF-15 into equine serum.
Sample material	Serum and Li heparin, K <sub>2</sub> EDTA and K <sub>3</sub> EDTA plasma
Sample stability	Stable for 8 days at 20 – 25 °C, 14 days at 2 – 8 °C, 12 months at -20 °C (±5 °C). Freeze only once.
Reagent stability	<b>cobas e 411, e 601, e 602, E170:</b> After opening at 2 – 8 °C: 12 weeks on the analyzers: 8 weeks; <b>cobas e 801:</b> on the analyzer: 16 weeks
Sample volume	<b>cobas e 411, e 601, e 602, E170:</b> 35 µL <b>cobas e 801:</b> 21 µL
LoB, LoD, LoQ* (Specification)	LoB: ≤ 350 pg/mL; LoD: ≤ 400 pg/mL; LoQ: ≤ 400 pg/mL
Measuring range	400 – 20,000 pg/mL
Dilution	1:5 (concentration of diluted sample must be >3500 pg/mL) using Diluent MultiAssay

\* LoB = Limit of Blank; LoD = Limit of Detection; LoQ = Limit of Quantitation (20 % CV)

## Order information

Elecsys® GDF-15	07125933 190 (E170, <b>e 411, e 601, e 602</b> ) 07028172 190 ( <b>e 801</b> )	100 tests per rackpack 100 tests per rackpack
CalSet GDF-15	07125941 190 (E170, <b>e 411, e 601, e 602, e 801</b> )	4 × 1.0 mL each of CalSet GDF-15 Level 1 and 2
PreciControl Cardiac II	04917049 190 (E170, <b>e 411, e 601, e 602, e 801</b> )	4 × 2.0 mL each of PreciControl Cardiac II 1 and 2
Diluent MultiAssay	03609987 190 (E170, <b>e 411, e 601, e 602</b> ) 07299010 190 ( <b>e 801</b> )	2 × 16 mL sample diluent 45.2 mL sample diluent

## Broad cardiac menu available on cobas e analyzers

CK-MB
Digitoxin
Digoxin
GDF-15
Myoglobin
NT-proBNP
Troponin I
Troponin T hs