



Advances in the management of suspected preeclampsia with Elecsys[®] sFlt-1/PIGF ratio





“Preeclampsia is a common and potentially serious condition that presents a continuing challenge to clinicians due to the variable features and lack of diagnostic tests.”

Prof. Andrew Shennan, King's College London, London, UK

Preeclampsia

A serious pregnancy complication

Introduction to preeclampsia¹

Preeclampsia is a serious complication in pregnancy which affects both the mother and the unborn child. Women with preeclampsia develop high blood pressure and high protein in their urine.

Preeclampsia is defined by the new onset of elevated blood pressure and proteinuria after 20 weeks of gestation. It is considered severe if blood pressure and proteinuria are increased substantially or symptoms of end-organ damage (including fetal growth restriction) occur. There is no single reliable, cost-effective screening test for preeclampsia, and there are no well-established measures for primary prevention. Management before the onset of labour includes close monitoring of maternal and fetal status.²

The majority of cases develop in healthy women bearing their first child. In addition several medical conditions are associated with an increased preeclampsia risk such as chronic hypertension, diabetes and renal disease. The cause of preeclampsia is not fully understood, but there is growing evidence that angiogenic factors such as placental growth factor (PlGF) and soluble fms-like tyrosine kinase-1 (sFlt-1) play a major role in the development of preeclampsia.

Preeclampsia may develop from 20 weeks of gestation until 48 hours after delivery. It is most commonly diagnosed after 32 weeks of gestation. Early-onset disease (20 – 33 + 6 weeks) is associated with particularly serious threats for the mother and fetus.

Complications

Preeclampsia is associated with increased mortality and morbidity. Women are often unaware that they have preeclampsia, even when it is life-threatening.

Amongst others, the most common complications of preeclampsia include:

For the mother:

- Eclampsia – preeclampsia in combination with generalized seizures
- Convulsions
- Kidney damage / Kidney failure
- Abruption placentae
- Antepartum hemorrhage
- Cerebrovascular bleeding

For the fetus:

- Fetal growth retardation
- Low birth weight
- Kidney damage / Kidney failure
- Premature birth
- Antepartum hemorrhage
- Stillbirth
- Respiratory distress syndrome
- Necrotizing enterocolitis




Incidence

Preeclampsia is the most common hypertensive disorder during pregnancy. It occurs in 3–5% of pregnancies and is defined by maternal hypertension with proteinuria.³

Preeclampsia has the greatest effect on maternal and infant outcome; it is a leading cause of preterm birth and consequent neonatal morbidity and mortality.

Hospital costs for both mother and baby⁴

- Preeclampsia without eclampsia
- Delivery with concomittant complications
- Normal delivery



Health Economics impact of Preeclampsia

Preeclampsia puts a major financial burden on the health care cost for pregnancies. In 2005 the average cost estimate for preeclampsia alone (excluding normal delivery costs) was GBP 9,009 per pregnancy.⁴

With an estimated 8.5 million women affected by preeclampsia every year, the estimated annual cost of preeclampsia worldwide is GBP 76.6 billion (based on the 2005 estimate). Costs might still be an underestimation due to the lack of proper diagnosis.⁵

Preeclampsia diagnosis and care is not optimal even in developed countries. Between 2006 and 2008 in the UK, 20 out of 22 deaths linked to preeclampsia involved substandard care. The substandard care in 63% of these deaths was categorized as major and they were described as “undoubtedly avoidable”.⁶



Pathogenesis


Hypertension and proteinuria are the diagnostic criteria for preeclampsia but they are only symptoms of the pathophysiologic changes that occur in the disorder. sFlt-1 and PlGF are indicators of the endothelial dysfunction associated with preeclampsia.^{7,8}

Aid in diagnosis

Diagnosis of preeclampsia is currently based on variable clinical parameters. Complications of the disease may be serious even when hypertension and/or proteinuria are mild.

Currently there is no single, objective laboratory test for the diagnosis of preeclampsia.² There is a need for a rapid and accurate aid in diagnosing this common and potentially serious condition to facilitate effective clinical management and to improve outcome for mother and fetus.^{8,9}

Preeclampsia landmarks



Each landmark demands a change in place and pace of care (frequency of screening)

Diagnosis of preeclampsia

Current standard of diagnosing preeclampsia²

Diagnosis of preeclampsia is not always easy. For the diagnosis of preeclampsia, both hypertension and proteinuria must be present.



Preeclampsia

- Blood pressure: 140 mmHg or higher systolic or 90 mmHg or higher diastolic after 20 weeks of gestation in a woman with previously normal blood pressure. Systolic increased > 30 mmHg or diastolic increased > 15 mmHg in a patient with preexisting chronic hypertension
- Proteinuria: 0.3 g or more of protein in a 24-hour urine collection

Severe preeclampsia

- Blood pressure: 160 mmHg or higher systolic or 110 mmHg or higher diastolic on two occasions at least six hours apart in a woman on bed rest
- Proteinuria: 5 g or more of protein in a 24-hour urine collection or 3+ or greater on urine dipstick testing of two random urine samples collected at least four hours apart
- Other features: oliguria (less than 500 mL of urine in 24 hours), cerebral or visual disturbances, pulmonary edema or cyanosis, epigastric or right upper quadrant pain, impaired liver function, thrombocytopenia, intrauterine growth restriction



NICE Antenatal Care Guideline 2008⁹

The NICE guidelines recommend that the following risk factors for preeclampsia should be sought at the first visit: age 40+ years, nulliparity, pregnancy interval of more than 10 years, family history of preeclampsia, BMI 30+ kg/m², pre-existing vascular disease such as hypertension, pre-existing renal disease and multiple pregnancy. According to NICE, more frequent blood pressure measurements should be considered for any women who have any of the above risk factors.



The preeclampsia community guideline (PRECOG)¹⁰

Refinements to the NICE guidelines, recommending that specialist referral is offered if **any risk factor** from the following list is present:

- **Previous preeclampsia**
- **Multiple pregnancy**
- **Long term medical condition**
 - **Hypertension**
 - **Renal disease**
 - **Diabetes**
 - **Antiphospholipid antibodies**

Refinements to the NICE guidelines, recommending that specialist referral is offered if **any two risk factors** from the following list are present:


- **First pregnancy**
- **≥ 10 years since the last baby**
- **Age ≥ 40 years**
- **BMI ≥ 35**
- **Family history of preeclampsia (mother or sister)**
- **Booking diastolic blood pressure ≥ 80 mmHg**
- **Proteinuria at booking ≥ 0.3 g / 24 hours**

Today's challenges in the clinical management of preeclampsia

The management of preeclampsia is troublesome

Clinical management of Preeclampsia can be challenging due to:

- A lack of effective treatment options other than delivery¹
- The necessity of a reliable tests for assessing disease severity and progression¹, and for predicting preeclampsia related complications¹¹
- The necessity to balance maternal and fetal risks¹², which may be in opposition
- The need for a reliable identification of high-risk preeclampsia patients to be referred to a specialized perinatal care center¹³



Diagnostic gold standard is not perfect

Current 'gold standard' for the diagnostic assessments include the determination of blood pressure and proteinuria.^{14,15,16}

However, **the current gold standard measurements have low sensitivity and specificity for predicting the disease and maternal and perinatal outcomes.**¹ Prediction and prognosis of preeclampsia remains an important unmet medical needs that is particularly critical given the unpredictable nature of this condition.

Also diagnosis of preeclampsia may be challenging when the clinical picture is not straightforward.



“Our understanding of the pathophysiology of preeclampsia, including the role of the placental factors sFlt-1 and PlGF, has improved. With this better understanding comes the opportunity to improve the way we diagnose this common and sometimes serious condition.”

Dr. Nadia Berkane, Hôpital Tenon, Paris, France

Elecsys® sFlt-1/PIGF


The first fully automated tests for preeclampsia

Improved understanding leads to diagnostic advances: sFlt-1 and PIGF


Preeclampsia may be caused by an imbalance of angiogenic factors. It has been demonstrated that high serum levels of sFlt-1, an anti-angiogenic protein, and low levels of PIGF, a pro-angiogenic protein, predict subsequent development of preeclampsia. In the absence of glomerular disease leading to proteinuria, sFlt-1 is too large a molecule to be filtered into the urine, while PIGF is readily filtered. The hypoxic placenta, which is commonly found in preeclampsia, produces sFlt-1.^{7,17,18,19,20}

These pathological changes lead to a vasospasm which is responsible for reduced perfusion of the placenta.


Hypoxic placenta



Normal pregnancy




Preeclampsia




In a multicenter case-control study the Elecsys sFlt-1 and Elecsys PIGF assays were tested in parallel on samples from 468 pregnant women with normal pregnancy outcome (no preeclampsia, no HELLP syndrome, no intrauterine growth restriction) and 234 patients with preeclampsia / HELLP syndrome. All pregnancies were singleton pregnancies.^{23,33} Different sets of cut-offs are suggested for early-onset and late-onset preeclampsia for use as an aid in the diagnosis of preeclampsia.

sFlt-1/PIGF ratio for aid in diagnosis of preeclampsia³³



PE/HELLP = preeclampsia / HELLP syndrome; sFlt-1/PIGF ratio cut-offs for early gestational phase: 33/85; sFlt-1/PIGF ratio cut-offs for late gestational phase: 33/110

sFlt-1/PIGF ratio for aid in diagnosis of preeclampsia³³



P < 0.001 multigroup comparison PE subgroups

“Preeclampsia by itself cannot be treated, but the clear stratification of risk can trigger concrete actions, such as the close monitoring of the mother and fetus as well as the referral to a specialist delivery unit offering intensive care.”

Prof. Holger Stepan, University of Leipzig, Germany

Improve the management of suspected preeclampsia patients

Patients with signs and/or symptoms of preeclampsia are sometimes difficult to be managed. Often, a suspicion of preeclampsia triggers escalation to a higher level of care. In fact, appropriate and timely referral to specialized centers reduces perinatal morbidity and mortality by 20%.¹³

The prognostic performance of the current diagnostic standard in determining which women will develop preeclampsia and how the disease will progress is quite poor.¹ As a consequence, many pregnant women with signs and/or symptoms of preeclampsia are often unnecessarily hospitalized for observation, resulting in stress for the expectant woman and significant additional costs to pregnancy care.

The measurement of the Elecsys[®] sFlt-1 / PIGF ratio is a reliable tool to identify the patients that are at high risk to develop preeclampsia requiring a closer monitoring. On the other hand, the Elecsys sFlt-1/PIGF test allows you to confidently send home patients for one week that are not going to develop the disease.^{22,23}

Improve outcome for mother & child through effective clinical management

In addition to supporting the prediction of preeclampsia Elecsys immunoassays help to optimise the clinical management. Following the diagnosis of preeclampsia an assessment is needed to grade the severity of the disease to determine whether conservative or active management is appropriate.

Decisions are needed as to whether urgent admission, hospital assessments or monitoring are appropriate.²

With the Elecsys sFlt-1/PIGF ratio the physician does not just have to rely on the degree of hypertension, the degree of proteinuria and the presence or absence of symptoms.

Addressing unmet medical needs

Medical value of Elecsys[®] preeclampsia immunoassays

Aid in short-term prediction of preeclampsia^{22,23}

A recent multicenter prospective study – PROGNOSIS (Prediction of short-term outcome in pregnant women with suspected preeclampsia study) – evaluated the use of the Elecsys sFlt-1/PIGF ratio in the short-term prediction of preeclampsia/ eclampsia/HELLP syndrome in pregnant women with suspected preeclampsia. The PROGNOSIS Study collected samples and clinical data from 1,273 pregnant women with clinical suspicion of preeclampsia between gestational weeks 24 + 0 days – 36 + 6 days at 30 study sites in different global regions.

One single cut-off of 38 for sFlt-1/PIGF ratio was identified in the PROGNOSIS Study:²³

- **sFlt-1/PIGF ratio < 38: Rule-out preeclampsia for 1 week**
- **sFlt-1/PIGF ratio ≥ 38: Rule-in preeclampsia within 4 weeks**

Short-term prediction of preeclampsia Rule-out within 1 week (n = 1,050)²³

sFlt-1/PIGF ratio	< 38
NPV (95% CI)	99.1% (98.2 – 99.6)
Sensitivity (95% CI)	85.7% (72.8 – 94.1)
Specificity (95% CI)	79.1% (76.5 – 81.6)

Short-term prediction of preeclampsia Rule-in within 4 weeks (n = 1,050)²³

sFlt-1/PIGF ratio	≥ 38
PPV (95% CI)	38.6% (32.6 – 45.0)
Sensitivity (95% CI)	70.3% (61.9 – 77.8)
Specificity (95% CI)	83.1% (80.5 – 85.5)

CI = Confidence Interval

NPV = Negative Predictive Value

PPV = Positive Predictive Value




Elecsys® sFlt-1/PIGF ratio as tool to support physician decision-making

When a pregnant woman presents with signs and symptoms of preeclampsia, the physicians often face a challenging decision making process. There are different components to be considered: the necessity to balance maternal and fetal risks¹², the need to best allocate healthcare resources to the patients at the highest risk but also the desire to avoid unnecessary stress for the patients.

The Preeclampsia Open Study (PreOS) assessed the influence of the Elecsys sFlt-1/PIGF test on the decision making of the physicians to hospitalize patients with suspicion of preeclampsia. PreOS collected data from 209 patients in a multicenter, prospective and non-interventional setting.^{24,25} The Elecsys sFlt-1/PIGF ratio was found to influence decision-making for hospitalisation in suspected preeclampsia and the changed decisions were in concordance with the incidence of major clinical outcomes (e.g. preeclampsia).

Preeclampsia outcome by hospital admission²⁵






Differential diagnosis of preeclampsia is often complicated

The clinical presentation of preeclampsia and subsequent clinical course of the disease can vary tremendously. The tools currently available to diagnose preeclampsia include measuring blood pressure and assessing proteinuria. However, these have low sensitivity and specificity in terms of assessing disease severity or predicting the course of the disease.¹

Underlying chronic disorders, such as renal diseases or autoimmune disorders can mimic the preeclampsia phenotype, making diagnosis difficult. Other hypertensive disorders in pregnancy include gestational hypertension and chronic hypertension. Pregnant women with chronic hypertension are at higher risk of developing superimposed preeclampsia than normotensive women.^{1,13,25}

Hypertensive pregnancy disorders: classification and diagnostic criteria¹⁹



* The diagnosis of eclampsia, a convulsive form of preeclampsia, is based on new-onset seizures, in the absence of a previous history of a seizure disorder. (GW = gestational weeks; HTN = hypertension)


“We hope that early identification of women at high risk of developing preeclampsia will allow healthcare professionals to prevent the most serious effects of the disease and avoid unnecessary expenditure by healthcare systems on excessive medical treatment or unnecessary hospital admission prompted by inadvertently positive diagnoses based on current standard of clinical practice”

Prof. Harald Zeisler, University of Vienna, Austria

Angiogenic factors can support in the differential diagnosis of preeclampsia

A study showed that the measurement of the sFlt-1/PIGF ratio can differentiate between different forms of hypertensive disorders. Women with preeclampsia or HELLP syndrome had significantly higher sFlt-1/PIGF ratios than women with gestational hypertension, chronic hypertension or no hypertensive disorder at all ($p < 0.001$).¹³

Another study showed that the sFlt-1/PIGF ratio may facilitate the diagnosis of superimposed preeclampsia in women with chronic hypertension.²⁶



sFlt-1/PIGF ratio in PE/HELLP, GH, chrHTN, and healthy controls¹³
 (PE = preeclampsia; HELLP = Hemolysis, Elevated Liver Enzymes, Low Platelets; GH = gestational hypertension; chrHTN = chronic hypertension)



Combination of angiogenic factors with Doppler sonography

Preeclampsia is characterised by an abnormal perfusion of the uterine arteries. Doppler sonography is often used as part of the clinical examination for patients with suspected preeclampsia, however, it has limited predictive value.²⁷

Combining Doppler sonography and biomarkers can improve the positive predictive value (PPV) for preeclampsia. Combination with biomarkers has been shown to improve the predictive performance of Doppler sonography.^{20,27}

For example, 2nd trimester sFlt-1 measurements improve the PPV for Doppler sonography from:

- 33 % to 50 % for all cases of Preeclampsia
- 31 % to 56 % for cases of Preeclampsia where delivery before 34 weeks was required²⁸

Use of angiogenic biomarkers in the clinical management of preeclampsia

Preeclampsia is associated with a number of serious adverse events, including:¹¹

- Maternal acute renal failure, liver dysfunction, seizures and cerebral accidents
- Fetal growth restriction
- Maternal or fetal mortality

Clinical criteria alone (blood pressure and proteinuria) may be inadequate to predict adverse outcomes.

Recent studies showed that a high sFlt-1/PlGF ratio and a more rapid elevation in the sFlt-1/PlGF ratio are associated with a significantly increased risk for an immediate delivery.^{11,13,20,29}

Health Economic impact of preeclampsia testing

Clinical use of the sFlt-1/PlGF ratio at triage could enable reductions in direct hospital costs and resource use

The use of the sFlt-1/PlGF assays could allow cost-savings


A correct short-term prediction and diagnosis of preeclampsia might allow a reduction of maternal and fetal morbidity as well as mortality. Timely referral of the patients at highest risk to a perinatal care centre reduces perinatal morbidity and mortality by 20%.¹³ In addition, quick and reliable detection of preeclampsia facilitates prompt intervention with steroids for fetal lung maturation, magnesium sulphate for seizure prophylaxis and antihypertensive therapy.

An improved prediction and diagnosis can allow a reduction of inappropriate discharges as well as a reduction of unnecessary hospitalizations, therefore a reduction of the health care burden.¹¹


In a recent model, adding the Elecsys® sFlt-1/PlGF ratio to the standard diagnostic method improved risk stratification with a significant reduction of false positive and false negative diagnosis. This could enable reductions in direct hospital costs and resource use with savings of 540 – 1,215 USD per patient.³⁰

By using the novel preeclampsia test in the UK, the NHS could save GBP 730 million annually and in Germany, national savings could reach EUR 436 million annually.^{31,32}


Standart approach³⁰



New approach³⁰
(adding sFlt-1/PlGF ratio)



Ideal case budget impact by implementation of sFlt-1/PlGF³⁰



Novel innovative Elecsys sFlt-1/PIGF ratio test

Precise, consistent, reliable

sFlt-1/PIGF ratio cut-offs^{23,24,33}

Early-onset preeclampsia – gestational week 20 – 33 + 6 days

sFlt-1/PIGF ≥ 85	▶ Diagnosis	▶ 99.5 % specificity the woman has preeclampsia Sensitivity: 88.0 %
sFlt-1/PIGF < 85 ≥ 38	▶ Prediction rule-in within next 4 weeks	▶ 38.6 % PPV the woman is at high risk to develop preeclampsia within the next 4 weeks
sFlt-1/PIGF < 38	▶ Prediction rule-out for the next 1 week	▶ 99.1 % NPV the woman will not develop preeclampsia for the next 1 week

Make sure you don't miss preeclampsia diagnosis

Roche Elecsys® sFlt-1/PIGF assay helps identifying preeclampsia patients at the highest need of intensified care in daily management and it provides reliable results in critical situations, when only a non-specific clinical picture is available. With the support of the Elecsys sFlt-1/PIGF assay, timely decisions regarding treatment and disease management are possible, thereby reducing danger to mother and fetus.

The sFlt-1/PIGF ratio allows gynaecologists to focus on the high-risk patients and reassure the patients not going to develop pre-eclampsia/eclampsia/HELLP syndrome for one week.^{23,24}

Late-onset preeclampsia – gestational week 34 to end of pregnancy

sFlt-1/PIGF ≥ 110	▶ Diagnosis	▶ 95.5 % specificity the woman has preeclampsia Sensitivity: 58.2 %
sFlt-1/PIGF < 110 ≥ 38	▶ Prediction rule-in within next 4 weeks	▶ 38.6 % PPV the woman is at high risk to develop preeclampsia within the next 4 weeks
sFlt-1/PIGF < 38	▶ Prediction rule-out for the next 1 week	▶ 99.1 % NPV the woman will not develop preeclampsia for the next 1 week

Peace of mind for physicians to timely and correctly diagnose preeclampsia

- a disease where the clinical picture is not always clear
- a disease where the onset is unpredictable
- a disease where the assessment of the severity/prognosis is difficult
- a disease where timely decision can be critical both for mother and fetus

References

- 1 Verlohren, S., Stepan, H., Dechend, R. (2012). Angiogenic growth factors in the diagnosis and prediction of pre-eclampsia. *Clin Sci*. 122(2), 43-52.
- 2 Magee, L.A., Pels, A., Helewa, M., Rey, E., von Dadelszen, P., on behalf of the Canadian Hypertensive Disorders of Pregnancy (HDP) Working Group (2014). Diagnosis, evaluation, and management of the hypertensive disorders of pregnancy. *Pregnancy Hypertension* 4(2), 105-145.
- 3 Chaiworapongsa, T., Chaemsaitong, P., Yeo, L., Romero, R. (2014). Pre-eclampsia part 1: current understanding of its pathophysiology. *Nat Rev Nephrol*. 10(8), 466-80.
- 4 Meads, C.A., Crossen, J.S., Meher, S., et al. (2008). Methods of prediction and prevention of pre-eclampsia: systematic reviews of accuracy and effectiveness literature with economic modelling. *Health Technol Assess*. 12(6), iii-iv.
- 5 Anderson, U.D., Olson, M.G., Kristensen, K.H., Åkerström, B., & Hansson, S.R. (2012) Review: Biochemical markers to predict preeclampsia. *Placenta*. 33(suppl), S42-7.
- 6 Shennan, A.H., Redman, C., Cooper C., & Milne, F. (2012). Are most maternal deaths from pre-eclampsia avoidable? *Lancet*. 379(9827), 1686-1687.
- 7 Andraweera, P.H., Dekker, G.A., Roberts C.T. (2012). The vascular endothelial growth factor family in adverse pregnancy outcomes. *Hum Reprod Update* 18(4), 436-57.
- 8 Rana, S., Karumanchi, S.A., Lindheimer, M.D. (2014). Angiogenic factors in diagnosis, management, and research in preeclampsia. *Hypertension* 63(2), 198-202.
- 9 National Institute for Health and Clinical Excellence (NICE) (2008). Antenatal care-routine care for the healthy pregnant woman. *NICE clinical guideline 62. 2nd edition. London: RCOG Press.*
- 10 Milne, F., Redman, C., Walker, J., Baker, P., Bradley, J., Cooper, C., et al. (2005). The pre-eclampsia community guideline (PRECOG): how to screen for and detect onset of pre-eclampsia in the community. *Br Med J*. 30, 576-580.
- 11 Rana, S., Powe, C.E., Salahuddin, S., Verlohren, S., Perschel, F.H., et al. (2012). Angiogenic factors and the risk of adverse outcomes in women with suspected preeclampsia. *Circulation*, 125(7), 911-919.
- 12 Chaiworapongsa, T., Chaemsaitong, P., Korzeniewski, S.J., Yeo, L., Romero R. (2014). Pre-eclampsia part 2: prediction, prevention and management. *Nat Rev Nephrol*. 10(9):531-40.
- 13 Verlohren, S., Herraiz, I., Lapaire, O., Schlembach, D., Moertl, et al. (2012). The sFlt-1/PlGF ratio in different types of hypertensive pregnancy disorders and its prognostic potential in preeclamptic patients. *Am J Obstet Gynecol*. 206, 58.e1-8.
- 14 National Institute for Health and Clinical Excellence (NICE) (2011). Hypertension in pregnancy. *NICE clinical guideline 107, August 2010, revised reprint January 2011, London: RCOG Press.*
- 15 ACOG Task Force on Hypertension in Pregnancy (2013). Hypertension in Pregnancy. *Obstet Gynecol*. 122, 1122-1131.
- 16 Williams, D., Craft, N. (2012). Pre-eclampsia. *Br Med J*. 345, e4437.
- 17 Chaiworapongsa, T., Chaemsaitong, P., Korzeniewski, S.J., Yeo, L., Romero, R. (2014). Pre-eclampsia part 2: prediction, prevention and management. *Nat Rev Nephrol*. 10(9), 531-540.
- 18 Levine, R.J., Maynard, S.E., Qian, C., Lim, K.H., England, L.J., et al. (2004). Circulating angiogenic factors and the risk of preeclampsia. *N Engl J Med*. 350, 672-683.
- 19 Garovic, V.D. (2012). The role of angiogenic factors in the prediction and diagnosis of preeclampsia superimposed on chronic hypertension. *Hypertension* 59, 555-557.
- 20 Gómez-Arriaga, P.I., Herraiz, I., López-Jiménez, E.A., Escribano, D., Denk, B., Galindo, A. (2014). Uterine artery Doppler and sFlt-1/PlGF ratio: prognostic value in early-onset pre-eclampsia. *Ultrasound Obstet Gynecol*. 43(5), 525-32.
- 21 Verlohren, S., Galindo, A., Schlembach, D., Zeisler, H., Herraiz, I., et al. (2010). An automated method for the determination of the sFlt-1/PlGF ratio in the assessment of preeclampsia. *Am J Obstet Gynecol*. 202(2), 161.e1-161.e11.
- 22 Hund, M., Allegranza, D., Schoedel, M., Dilba, P., Verhagen-Kamerbeek, W. & Stepan, H., (2014) Multicenter prospective clinical trial to evaluate the prediction of short-term outcomes in pregnant women with suspected preeclampsia (PROGNOSIS): study protocol. *BMC Pregnancy and Childbirth*. 14, 324.
- 23 Elecsys® sFlt-1 package insert (Nov 2014) and Elecsys® PlGF package insert (Nov 2014). *Roche Diagnostics Documentation, Basel.*
- 24 Hund, M., Verhagen-Kamerbeek, W.D.J., Reim, M., Messinger, D., van der Does, R. & Stepan, H. (2014). Influence of the sFlt-1/PlGF ratio on clinical decision-making in women with suspected preeclampsia – the PreOS study protocol. *Hypertension in Pregnancy*. Accepted for publication.
- 25 Klein, E., Schlembach, D., Ramoni, A., Langer, E., Bahlmann, F., van der Does, R., Messinger, D., Verhagen-Kamerbeek, W.D.J., Reim, M., Hund, M. & Stepan, H. (2014). PreOS (Preeclampsia Open Study), a multicenter, prospective, non-interventional study evaluating the influence of the sFlt-1/PlGF ratio on physician decision-making in pregnant women with suspicion of preeclampsia. *20th World Congress on Contraversies in Obstetrics, Gynecology & Infertility (COGI).*
- 26 Perni, U., Sison, C., Sharma, V., Helseth, G., Hawfield, A. et al. (2012). Angiogenic factors in superimposed preeclampsia: a longitudinal study of women with chronic hypertension during pregnancy. *Hypertension* 59(3), 740-746.
- 27 Giguère, Y., Charland, M., Bujold, E., Bernard, N., Grenier, S., et al. (2010). Combining biochemical and ultrasonographic markers in predicting preeclampsia: a systematic review. *Clin Chem*. 56(3), 361-375.
- 28 Stepan, H., Unversucht, A., Wessel, N., & Faber, R. (2007). Predictive value of maternal angiogenic factors in second trimester pregnancies with abnormal uterine perfusion. *Hypertension*. 49(4), 818-824.
- 29 Álvarez-Fernández, I., Prieto, B., Rodríguez, V., Ruano, Y., Escudero, A.I., Álvarez, F.V. (2014). New biomarkers in diagnosis of early onset preeclampsia and imminent delivery prognosis. *Clin Chem Lab Med*. 52(8):1159-1168.
- 30 Schnettler, W.T., Dukhovny, D., Wenger, J., Salahuddin, S., Ralston, S.J., & Rana S. (2013). Cost and resource implications with serum angiogenic factor estimation in the triage of pre-eclampsia. *BJOG*. 120(10), 1224-1232.
- 31 Hadker, N., Garg, S., Costanzo, C., Miller, J.D., Foster, T., et al. (2010). Financial impact of a novel pre-eclampsia diagnostic test versus standard practice: a decision-analytic modeling analysis from a UK healthcare payer perspective. *J Med Econ*. 13(4), 728-737.
- 32 Hadker, N., Garg, S., Costanzo, C., van der Helm, W., & Creeden, J. (2013) Are there financial savings associated with supplementing current diagnostic practice for preeclampsia with a novel test? Learnings from a modeling analysis from a German payer perspective. *Hypertens Pregnancy*. 32(2), 105-119.
- 33 Verlohren, S., Herraiz, I., Lapaire, O., Schlembach, D., Zeisler, H., et al. (2014) New gestational phase-specific cutoff values for the use of the soluble fms-like tyrosine kinase-1/placental growth factor ratio as a diagnostic test for preeclampsia. *Hypertension*. 63(2), 346-352.

COBAS, ELECSYS and LIFE NEEDS ANSWERS are trademarks of Roche.

©2014 Roche

Roche Diagnostics International Ltd
CH-6343 Rotkreuz
Switzerland
www.cobas.com